



Recombinant Escherichia coli O6:K15:H31 Anaerobic nitric oxide reductase transcription regulator norR (norR)

Product Code	CSB-MP606994EGY
Abbreviation	norR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TEH1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSFSVDVLAN IAIELQRGIG HQDRFQRLIT TLRQVLECD A SALLRYDSRQ FIPLAIDGLA KDVLGRRFAL EGHPRLEAIA RAGDVVRFPA DSELPDPYDG LIPGQESLKV HACVGLPLFA GQNLIGALT DGMQPDQFDV FSDEELRLIA ALAAGALSNA LLIEQLESQN MLPGDAAPFE AVKQTQMIGL SPGMTQLKKE IEIVAASDLN VLISGETGTG KELVAKAIHE ASPRAVNPLV YLNCAALPES VAESELF GHV KGAF TGAIN RSGKFEMADN GTLFLDEIGE LSLALQAKLL RVLQYGD IQR VGDDRSLRVD VRVLAATNRD LREEVLAGRF RADLFHRLSV FPLSVPLRE RGDDVILLAG YFCEQCRLRL GLSRVLSAG ARNLLQHYNF PGNVRELEHA IHRAVVLARA TRSGDEVILE AQHFAPPEVT LPPPEAAAVP VVKQNLREAT EAFQRETIRQ ALAQNHNNWA ACARMLET DV ANLHRLAKRL GLKD
Source	Mammalian cell
Target Names	norR
Protein Names	Recommended name: Anaerobic nitric oxide reductase transcription regulator norR
Expression Region	1-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.