



Recombinant Escherichia coli O6:K15:H31 Anaerobic glycerol-3-phosphate dehydrogenase subunit B (glpB)

Product Code	CSB-EP607006EGY
Abbreviation	glpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TFK1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p> MRFDTVIMGG GLAGLLCGLQ LQKHGLRCAI VTRGQSALHF SSGSLDLLSH LPDGQPVTEI HSGLESLRQQ APAHPYTLTG PQRVLDLACQ AQALIAESGA QLQGSVELAH QRITPLGTLR STWLSSPEVP VWPLPAKKIC VVGISGLMDF QAHLAAASLR ELDLKVETAE IELPELDVLR NNATEFRAVN IARFLDNEEN WPLLLDALIP VANTCEMILM PACFGLANDK LWHHWLNEKLP CSLMLLPTLP PSVLGIRLQN QLQRQFVRQG GVVWMPGDEVK KVTCKNGVVN EIWTRNHADI PLRPRFAVLA SGSFFSGGLV AERDGI REPI LGLDVLQTAT RGEWYKGDFF APQPWQQFGV TTDEALRPSQ AGQTIENTLFA IGSVLGGFDP IAQGCGGGVC AVSALHAAQQ IAQRAGGQQ </p>
Source	E.coli
Target Names	glpB
Protein Names	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit B Short name= Anaerobic G-3-P dehydrogenase subunit B Short name= Anaerobic G3Pdhase B EC= 1.1.5.3
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.