



Recombinant Escherichia coli O6:K15:H31

Probable endonuclease 4 (nfo)

Product Code	CSB-MP607008EGY
Abbreviation	nfo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TFT1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKYIGAHVSA AGGLANA AIR AAEIDATAFA LFTKNQRQWR AAPLTTQTID EFKAACEKYH YTSAQILPHD SYLINLGHPV AEALEKSRDA FIDEMQRCEQ LGLSLLNFHP GSHLMQISEE DCLARIAESI NIALDKTQGV TAVIENTAGQ GSNLGFKFEH LAAIDGVED KSRVGVCIDT CHAFAAGYDL RTPAECEKTF ADFARIVGFK YLRGMHLNDA KSTFGSRVDR HHSLGEGNIG HDAFRWIMQD DRFDGIPLIL ETINPDIWAE EIAWLKAQQT EKAVA
Source	Mammalian cell
Target Names	nfo
Protein Names	Recommended name: Probable endonuclease 4 EC= 3.1.21.2 Alternative name(s): Endodeoxyribonuclease IV Endonuclease IV
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.