



Recombinant Escherichia coli O6:K15:H31 tRNA 2-thiocytidine biosynthesis protein TtcA (ttcA)

Product Code	CSB-EP607021EGY-B
Abbreviation	ttcA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TI22
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQENQQITKK EQYNLNKLQK RLRRNVGEAI ADFNMIEEGD RIMVCLSGGK DSYTMLEILR NLQQSAPINF SLVAVNLDQK QPGFPEHVLP EYLETLGVEY KIVEENTYGI VKEKIPEGKT TCSLCSRLRR GILYRTATEL GATKIALGHH RDDILQTLFL NMFYGGKMKG MPPKLMSDDG KHIVIRPLAY CREKDIQRFA DAKAFPIIPC NLCGSQPNLQ RQVIADMLRD WDKRYPGRIE TMFSAMQNVV PSHLCDTNLF DFKGITHGSE VVNGGDLAFD REEIPLQPAG WQPEEDENQL DELRLNVVEV K</pre>
Source	E.coli
Target Names	ttcA
Protein Names	Recommended name: tRNA 2-thiocytidine biosynthesis protein TtcA
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.