



Recombinant Escherichia coli O6:K15:H31 Flagellar P-ring protein (flgI)

Product Code	CSB-YP607027EGY
Abbreviation	flgI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TIZ4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	E RIRDLTSVQG VRQNSLIGYG LVVGLDGTGD QTTQTPFTTQ TLNNMLSQLG ITVPTGTNMQ LKNVAVMVT ASLPPFGRQG QTIDVVVSSM GNAKSLRGGT LLMTPLKGVD SQVYALAQGN ILVGGAGASA GGSSVQVNQL NGGRITNGAV IERELPSQFG VGNTLNLQLN DEDFSMAQQI ADTINRVRGY GSATALDART IQVRVPSGNS SQVRFLADIQ NMQVNVTPQD AKVVINSRTG SVVMNREVTL DSCAVAQGNL SVTVNRQANV SQPDPFPGGG QTVVTPQTQI DLRQSGGSLQ SVRSSASLNN VVRALNALGA TPMDLMSILQ SMQSAGCLRA KLEII
Source	Yeast
Target Names	flgI
Protein Names	Recommended name: Flagellar P-ring protein Alternative name(s): Basal body P-ring protein
Expression Region	20-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.