



# Recombinant Escherichia coli O6:K15:H31 Ribosomal protein S12 methylthiotransferase RimO (rimO)

<b>Product Code</b>	CSB-MP607034EGY
<b>Abbreviation</b>	rimO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TJL3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSKVTPQPKI GFVSLGCPKN LVDSEILTE LRTEGYDVVP SYDDADMVIV NTCGFIDSAV QESLEAIGEA LNENGVIVT GCLGAKEDQI REVHPKVLEI TGPHSYEQVL EHVHHYVPKP KHNPFLSLVP EQGVKLTPRH YAYLKISEGC NHRCTFCIIP SMRGDLVSRP IGEVLSEAKR LVDAGVKEIL VISQDTSAYG VDVKHRTGFH NGEPVKTSMV SLCEQLSKLG IWTRLHYVYP YPHVDDVIPL MAEGKILPYL DIPLQHASPR ILKLMKRPGS VDRQLARIKQ WRKICPELTL RSTFIVGFPG ETEEDFQMLL DFLKEARLDR VGCFKYSPVE GADANALPDQ VPEEVKEERW NRFMQLQQQI SAERLQEKVG REILVIIDEV DEEGAIGRSM ADAPEIDGAV YLNGETNVKP GDILRVKVEH ADEYDLWGSR V
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rimO
<b>Protein Names</b>	Recommended name: Ribosomal protein S12 methylthiotransferase RimO Short name= S12 MTTase Short name= S12 methylthiotransferase EC= 2.-.- Alternative name(s): Ribosome maturation factor RimO
<b>Expression Region</b>	1-441
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.