



Recombinant *Xenopus tropicalis* Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial (ecsit)

Product Code	CSB-BP607141XBF
Abbreviation	ecsit
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0V9C9
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	TGN NAQKQVTCCL HTSPPPGPPQ ATSSAGSVTP YEDVFQREQR NKASFVEVLD LYCNRDVRRR GHVEIEAAL RWMTEFGVEK DLEVYNKLLD VFPKEVFVPQ NYIQRMFNHY PRQQECAIRV LEQMENYGIT PNKQTCFLLL QIFGARSHPI RKYQRMMYWF PRFKHTNPFV VPAELPADPV ELSRMCLHRI ANDLDSRVAV YQMPYTEKCD DGKEREHLHV VGIQSPSQVS LLAHNPSTP VYVEGPFPLW LRKTCVYYYI LRADPDPMRK EEEIDNERSL YYPLSLDLEL DRDLGDDDSF DVDDVEEGPV FAMCMAGSGD ELLLGKWIRG LQETNPILGQ TPVVFRLNSG PQELSVPQEK EAEPKGELED DEGEREPERR RMEQ
Source	Baculovirus
Target Names	ecsit
Protein Names	Recommended name: Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial
Expression Region	28-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.