



Recombinant Human Secernin-3 (SCRN3)

Product Code	CSB-MP607167HU
Abbreviation	SCRN3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q0VDG4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MEPFSCDTFV ALPPATVDNR IIFGKNSDRL YDEVQEVVYF PAVVHDNLGE RLKCTYIEID QVPETYAVVL SRPAWLWGAE MGANEHGVCI GNEAVWGREE VCDEEALLGM DLVRLGLERA DTAEKALNVI VDLLEKYGGQ GNCTEGRMVF SYHNSFLIAD RNEAWILETA GKYWAAEKVQ EGVRNISNQL SITTKIAREH PDMRNYAKRK GWWDGKKEFD FAAAYSYLDT AKMMTSSGRY CEGYKLLNKH KGNITFETMM EILRDKPSGI NMEGEFLTTA SMVSILPQDS SLPCIHFFTG TDPERSVFK PFIFVPHISQ LLDTSSPTFE LEDLVKKKSH FKPDRRHPLY QKHQQALEVV NNNEEKAKIM LDNMRKLEKE LFREMESILQ NKHLDVEKIV NLFPQCTKDE IQIYQSNLSV KVSS</p>
Source	Mammalian cell
Target Names	SCRN3
Protein Names	Recommended name: Secernin-3
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>