



Recombinant Human Secernin-3 (SCRN3)

Product Code	CSB-YP607167HU
Abbreviation	SCRN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0VDG4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPFSCDTFV ALPPATVDNR IIFGKNSDRL YDEVQEVVYF PAVVHDNLGE RLKCTYIEID QVPETYAVVL SRPAWLWGAE MGANEHGVC I GNEAVWGREE VCDEEALLGM DLVRLGLERA DTAEKALNVI VDLLEKYGGQ GNCTEGRMVF SYHNSFLIAD RNEAWILETA GKYWAAEKVQ EGVRNISNQL SITTKIAREH PDMRNYAKRK GWWDGKKEFD FAAAYSYLDT AKMMTSSGRY CEGYKLLNKH KGNITFETMM EILRDKPSGI NMEGEFLT TA SMVSILPQDS SLPCIHFFTG TPDERSVFK PFIFVPHISQ LLDTSSPTFE LEDLVKKKSH FKPDRRHPLY QKHQQALEVV NNNEEKAKIM LDNMRKLEKE LFREMESILQ NKHLDVEKIV NLFPQCTKDE IQIYQSNLSV KVSS
Source	Yeast
Target Names	SCRN3
Protein Names	Recommended name: Secernin-3
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.