



# Recombinant Arabidopsis thaliana UDP-glycosyltransferase 75C1 (UGT75C1)

<b>Product Code</b>	CSB-EP607299DOA
<b>Abbreviation</b>	UGT75C1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0WW21
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATSVNGSHR RPHYLLVTFP AQGHINPALQ LANRLIHHGA TVTYSTAVSA HRRMGPEPST KGLSFAWFTD GFDDGLKSFE DQKIYMSELK RCGSNALRDI IKANLDATTE TEPITGVIYS VLVPWVSTVA REFHLPTLL WIEPATVLDI YYYYFNTSYK HLFDVEPIKL PKLPLITTGD LPSFLQPSKA LPSALVTLRE HIEALETESN PKILVNTFSA LEHDALTSVE KLKMIPIGPL VSSSEGKTDL FKSSDEDYTK WLDSKLETSV IYISLGTHAD DLPEKHMEAL THGVLATNRP FLWIVREKNP EEKKNRFLE LIRGSDRGLV VGWCSQTAVL AHCAVGCFVT HCGWNSTLES LESGVPVVAF PQFADQCTTA KLVEDTWRIG VKVKVGEEGD VDGEEIRRCL EKVMMSGEEA EEMRENAEKW KAMAVIDAAAE GGPSDLNLKG FVDEDE
<b>Source</b>	E.coli
<b>Target Names</b>	UGT75C1
<b>Protein Names</b>	Recommended name: UDP-glycosyltransferase 75C1 EC= 2.4.1.- Alternative name(s): Anthocyanin 5-O-glucosyltransferase UDP glucose:anthocyanin 5-O-glucosyltransferase
<b>Expression Region</b>	1-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.