



Recombinant Arabidopsis thaliana Rhodanese-like domain-containing protein 11, chloroplastic (STR11)

Product Code	CSB-EP607300DOA-B
Abbreviation	STR11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0WWT7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MQAV DEDIDLKQMR DIAAAKKRWD GLLREGKVKL LTPREAGYAI SLSNKPLLDV RPSSERNKAW IKGSTWVPIF DNDDNLDAGT LSKKVTSFAM GGWWSGAPTL SFNRLFLSKV EEKFPKDSEL IVACQKGLRS LAACELLYNA GYENLFWVQG GLESAQDEDL VTEGVQPLKL AGIGGFSEFL GWTDQQRAQA AKEGWGYRLV YTARLFGVVL AADALFVGAQ QLGHYIQELR GH
Source	E.coli
Target Names	STR11
Protein Names	Recommended name: Rhodanese-like domain-containing protein 11, chloroplastic Alternative name(s): Sulfurtransferase 11 Short name= AtStr11
Expression Region	57-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.