



Recombinant Pre-mRNA-processing factor 19 homolog (T10F2.4)

Product Code	CSB-BP607335CXY
Abbreviation	T10F2.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10051
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSFVCGISGE LTEDPVVSQV SGHIFDRRLI VKFIAENGTD PISHGELSE QLVSLKSGGT GSAPRNVSGT SIPSLKMLQ DEWDTVMLNS FSLRQQQLQIA RQELSHSLYQ HDAACRVISR LSKELTAARE ALSTLKPHTS AKVDDDVSID ESEDQQGLSE AILAKLEEKs KSLTAERKQR GKNLPEGLAK TEELAELKQT ASHTGIHSTG TPGITALDIK GNLSLTGGID KTVVLYDYEK EQVMQTFKGH NKKINAVVLH PDNITAIAS ADSHIRVWSA TDSSSKAIID VHQAPVTDIS LNASGDYILS ASDDSYWAFS DIRSGKSLCK VSVEPGSQIA VHSIEFHPDG LIFGTGAADA VVKIWDLKNQ TVAAAFPGHT AAVRSIAFSE NGYYLATGSE DGEVKLWDLR KLKNLKTFAN EEKQPINSLS FDMTGTFLGI GGQKVQVLHV KSWSEVVSLs DHSGPVTGVR FGENARSLVT CSLDKSLRVF SF
Source	Baculovirus
Target Names	prp-19
Protein Names	Recommended name: Pre-mRNA-processing factor 19 homolog Alternative name(s): PRP19/PSO4 homolog
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.