



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C11D3.10 (SPAC11D3.10)

Product Code	CSB-MP607342SXV
Abbreviation	SPAC11D3.10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10089
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNALAKIKAF RSYVPLLQGN TKTIYLNQSF QAPMNMLVST ALQGYINEGL YNPHPKPMWK ERTEETRSL AKLLNASTKD SITFTRDTTE GLNLFQRSLK WKPGDNVIL DNEHPNQGF G WIALQNDGLE VRLVPNEGQY HANASTFPAY VDSRTKAIGL SSVMFHSGQK NDVKDIANAF RPKGIHVLAD LTQQVGLSKI DVQDLNVSAC AFSCHKGLGC PTGLGVLYVS PLAISELRST PPFVGGGAVE DFKEDLKLKL NAKYHQ SALR YEHTNNAYML ITALRAYLKF LLKVGISNVE RYLQGLGKDL IKELESLNVS VIGYKDFDKH SSHSYVLKIL NPEWFDLFRQ QGVCVSRFES GIRVSFGLYN TSKDIIKFIS VIRKGLALNI PLNIRPPQRI AVMDNPLVGI SDEAIAAMKL PANTLSNRPL AANI
Source	Mammalian cell
Target Names	SPAC11D3.10
Protein Names	Recommended name: Uncharacterized protein C11D3.10
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.