



Recombinant Schizosaccharomyces pombe

Uncharacterized protein C18G6.12c

(SPAC18G6.12c)

Product Code	CSB-EP607345SXV
Abbreviation	SPAC18G6.12c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10111
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAFSFQVSKM NSFDNYNAAI HPKHDLISVN APVFGKDAKL SPHLAIVLAD KGFMLEQVCL PWRYFTDRGF VVEFVIAESI FPPRPPKPND SYMTGWKSHV LGASKEILDI YSVLCTLNEF NNPCKYRSNG FTFDNFSAVF ITGGRNPFVR EMLEDPLLHA SLVPYIHSCR TIQPDEEVKD NRKIKVLGAI SQGAAAIYLA EPNINMKTTT IPAWMERSNS FLNPTPDNST YPYAATIIPK DKYVSGPYRR VAFTYQDENY YYISGRSNKD IGELSKKMYL MYKQAYKDLN ASLRERRRRS TSNRRNSGI
Source	E.coli
Target Names	SPAC18G6.12c
Protein Names	Recommended name: Uncharacterized protein C18G6.12c
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.