



Recombinant Schizosaccharomyces pombe Ribosomal RNA-processing protein 8 (rrp8)

Product Code	CSB-EP607376SXV-B
Abbreviation	rrp8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10257
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFKVDWDLGP VSKSASDQTK ESRKEKKRKK GERKNVGDKG EKLNEKVLKK AKSVTTNNSL KSEIKKEKSV PSIKEKNKGD AKHTKLTSLQ QKMKDKLDGA NFRWINEQLY TTESDKAVQM FKENPDLFDI YHAGFRYQVE GWPENPVDIF IQHLKIRFEH SNAKKKNNIV IADLGCGEAK IASTFRKSRS LQVHSFDLVA PNEHVVACDI ANVPMADETV DIAVFCLSLM GTNWQSFLKE AYRILKVGGL LWVAEIKSRF SDKSGEVFAK ELPKLGFKETK SIQLQNKMFT LFEFKKVPVH GKCEELPPIL SACIYKRR
Source	E.coli
Target Names	rrp8
Protein Names	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.-
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.