



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1F12.08 (SPAC1F12.08)

Product Code	CSB-EP607395SXV
Abbreviation	SPAC1F12.08
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10350
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGKTLIFIKL SSLNECIKYS TIISKIAREL SIHKGNELYI WVSCSEMCQK RTKHLFRELQ NALCELYLSS ATTNDSSLLS VVPTIVLFY WSNIQLGSIT DSWNSVYYSE NAKDQILNIP AKSFQSIGPD VLNEIHNIEH LPKKDESEDS GKEINNMVSA VGGTFDHLHV GHKVLLTLTA WFGVKEVIVG VSGDELLKKK VQKEFLENIQ KRKEEVSNFL HSIKEDINCR VVTIHDPFGP TITDAEIDSL IVSEETKTGA TAVQEERVKR GLPELSIYCI DLLYPAQELN LKDNNSLKVS STAIREEELAK LAYKNK
Source	E.coli
Target Names	SPAC1F12.08
Protein Names	Recommended name: Uncharacterized protein C1F12.08
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.