



Recombinant Protein mom-2 (mom-2)

Product Code	CSB-EP607409CXY-B
Abbreviation	mom-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10459
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	WWLLSK TDTSSANS GS SPILCKNVPG LTPQQKRMCH ENPNIKYLI SGLRSALHTC EYTFQREAWN CTLTLPGVGT SPLQIASRES AYYVAISAAG VSHSLARACS KGLIDDCGCG ETPQGGSGSVA VSQASSRSSS DFVWAGCSDN VKFGNTFGRK FVDQYDRQHA TEPRSQMNLH NNRVGRLLV NAMNKECKCH GVSGSCVTKT CWKVMPKFDE FASRLHQKYQ LAKLVTNNDQ KLTVRSSPSA GSSGRSERFA RNMDASSQM RNELIYLDAS PNYCAIDVKD RECGENCPNI CCGRGWRTRR EIVDEPCHCQ FVWCCEVKCK TCKKLVERNY CL
Source	E.coli
Target Names	mom-2
Protein Names	Recommended name: Protein mom-2
Expression Region	25-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.