



Recombinant Transcription factor elt-2 (elt-2)

Product Code	CSB-MP607427CXY
Abbreviation	elt-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10655
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDNNYNDNVN GWAEMEPSQP MGGLRLPTQN MDPPEQNNES QLSELPRMKI DNDYASPIER QSVITSGTNN YEPKVETVTS FFHTGIDYSN FGMLDQTTMQ PFYPLYSGIP VNTLGTFSGY TNSIYDKPSL YDPSIPTINI PSTYPTVAPT YECVKCSQSC GAGMKAVNGG MMCVNCSTPK TTYSPPVAYS TSLGQPPILE IPSEQPTAKI AKQSSKKSSS SNRGSNGSAS RRQGLVCSNC NGTNTTLWRR NAEGDPVCNA CGLYFKLHHI PRPTSMKKEG ALQTRKRKSK SGDSSTPSTS RARERKFERA SSSTEKAQRS SNRRAGSAKA DRELSTAAVA AATATYVSHA DLYPVSSAAV TLPDQTYSNY YQWNTAATAG LMMVPNDQNY VYAATNYQTG LRPADNIQVH VMPVQDDETK AAARDLEAVD GDS
Source	Mammalian cell
Target Names	elt-2
Protein Names	Recommended name: Transcription factor elt-2
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.