



Recombinant Protein dom-3 (dom-3)

Product Code	CSB-YP607428CXY
Abbreviation	dom-3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q10660
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MSHYGGNPRG NSSHQFGRKD FQQSDSKHIP KITGQPLPNE VQIPFDNMIY ETRNPPKFEK QAKFISEYCI NYDRKLQLGR MRAKKFHDKL PPNNLSLDLG KGFETFDPKE GDEKIIMLLE WINQMAPEGG RLKKVIHEAD VVCWRGLITK ICSTIYNKEE GWRIDAIKRK GVIFLCEEKT EQAKNRDLNQ TLDKRMSYW GHKFEQYVTL DDDAEQPPTS SPVTTKEEYA VVFRNDLQTD QRLNPEKRSI GIFYSGEVDC LDRHGNMLEL KTQKGELGHG FYKYSKSFKW WLQSSLVNVD HIIVGLRTQE GHVKSLTSLR TREIPQRASW NFRAGFEFLS TIFTYILNCL EKDGDACVIE YRSEMGIHKG IQMRRVPVDE FDFVPNEFLE KFC</p>
Source	Yeast
Target Names	dom-3
Protein Names	Recommended name: Protein dom-3
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>