



Recombinant *Loxodonta africana* F-actin-capping protein subunit alpha-2 (CAPZA2)

| | |
|--------------------------|---|
| Product Code | CSB-YP607446LJC |
| Abbreviation | CAPZA2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q108U5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Loxodonta africana</i> (African elephant) |
| Purity | >85% (SDS-PAGE) |
| Sequence | ADLEEQLSD EEKVRIAAKF ITHAPPGEFN EVFNDVRLLL NNDNLLREGA AHAFAYNLD QFTPVKIEGY EDQVLITEHG DLGNGKFLDP KNRICFKFDH LRKEATDPRP YEAENAVESW RTSVETALRA YVKEHYPNGV CTVYGKKIDG QQTIACIES HQFQAKNFWN GRWRSEWKFT ITPSTTQVVG ILKIQVHYE DGNVQLVSHK DIQDSLTVSN EVQTAKEFIK IVEAAENEYQ TAISENYQTM SDTTFKALRR QLPVTRTKID WNKILSYKIG KEMQNA |
| Source | Yeast |
| Target Names | CAPZA2 |
| Protein Names | Recommended name: F-actin-capping protein subunit alpha-2 Alternative name(s): CapZ alpha-2 |
| Expression Region | 2-286 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |