



Recombinant Nuclear hormone receptor family member nhr-2 (nhr-2)

Product Code	CSB-EP607447CXY-B
Abbreviation	nhr-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10902
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVPMNPELMH YPNQPHMMPG YPMINAYQQS YLQNPAQWTP EMLLCMQMQM QQMVPGTATN STYMPPQSKN TVTPNATSPL SQISDRDSGN DTISPPLTSQ NSATTPNYHN QPQMQLTPLQ QAKPEKEEIS PYAINNILSS TQSSPDNGNE ADDKENS SVR RNTFPASSAS TTKPRGTPYD RTNTPMTVVP QMPQFAPDFS NYNLTMNAWS EFFKKDRCMV CGDNSTGYHY GVQSCEGCKG FFRRSVHKNI AYVCTKGENC TFSYENCAAN RGRVTRCQAC RFAKCLAVGM NRDNVRVNKE TDKDVKPSVA SPNFEMTSQV KELTAAFVAN MPCSTHLTSG THAIGAIAKKF IESVPALSSL LPKDEKALEM SIQKVMGIL AIRAAFTFDP ITFYSCENPV NLLRGGIRNT VFNDCEVALL SGIHILQIVN GGVAEIFTSY CQGLRHQLSQ THLQEIGLCD RLLMRLGPYL NQ
Source	E.coli
Target Names	nhr-2
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-2
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.