



Recombinant *Oryza sativa* subsp. *japonica* Putative cyclin-D2-3 (CYCD2-3)

Product Code	CSB-EP607465OFG
Abbreviation	CYCD2-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10K98
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGFLCRRRSS AAPAESFLL CEEDSESVFG SDDDGVEETA TMAPELGKMM SLGFSASHHL GDGGGGGEEL VGSFMEKEVE QMVETARGEY LTKLSNGGIE LSCRIAIDW ICKVQAYYSF GPLCAYLAVN YLDRFLSSVE FSVTNDMPWM QQLLIVACLS LAAKMEETAA PGTLDLQVCN PEYVFDAETI HRMEIIVLTT LKWRMQAVTP FTYIGHFLDK INEGNRITSE LISRCTEIL STMKATVFLR FRPSEIATAV ALSVVADGGR VLDFGGVLES SKLPVDKDNV GRCHQAMQEM ALVMQNSTAS PSDSLCYDMT NPVLDCHILR RRECVSLQSG PRDELEMGWA PFAITVSEKE ESERRCDDGW RRGGEAAQYC EQADEQRRST DGRAGVRLEE QAGGE
Source	E.coli
Target Names	CYCD2-3
Protein Names	Recommended name: Putative cyclin-D2-3 Alternative name(s): G1/S-specific cyclin-D2-3 Short name= CycD2;3
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.