



Recombinant *Oryza sativa* subsp. *japonica* Putative cyclin-F3-1 (CYCF3-1)

Product Code	CSB-YP607479OFG
Abbreviation	CYCF3-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10Q63
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEAAAAAAAAE EEAGNPDGVE GAAVAVAPE AAAEGPSEPN AGEASREPDA GQASREPGAA GPSREPDVAG PSREPDAAGP SREPGAAGGS REPGAAGGSR QPVPDAAQLA VVPYVEDIDR YLRSLEAEQT RRPMINYVQE IQGGIINMDV RGILVDWMAD VAYVFNLQEE TLHHAVSYVD RFLSKIAFPG DKLKLLGTTA LfvaskyEEI HPPHVRNFSA VTVNTYTTQQ VSKMELDILR FLNFDVGSPT VITFLRKFLT SCCGGNNSN RKLELMCNYL AELSLDDYY IRFLPSIVAA ACLFVGKFTL NPNTRPWFSG VSTITPPENI KGGVEKYMVS RIYMCVFDLP MLFLMETWSS SGVSNHTKLQ LKQNMMELEPI ENPTFISAS
Source	Yeast
Target Names	CYCF3-1
Protein Names	Recommended name: Putative cyclin-F3-1 Short name= CycF3;1
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.