



Recombinant Uncharacterized protein C03F11.2 (C03F11.2)

Product Code	CSB-MP607535CXY
Abbreviation	C03F11.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11123
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MNLTASVSMV YDGGHYIVVD SPSATDVHSK ELMLKGIASR NIVPGEIQYV VTTHGHPDHF GQGNFFPNAR HFFGSYEYSD TNFISTELHT KDIMQLTKNV QLWNTPGHTA QDVTVMVHNV SCCGIIAVAG DLFYNEEDAN EAAGIWFQEA WNPIIGKISR NKVICYADYV IPGHGKLFRI TQEMKNGADC FTKYETTLDD NKTQVQSSPQF QNSIEQVIRK PAETLFSSAN SVSFNYIGSP STSASSSTRA IETTTAYSIDS TTTTFYTVDL TAETPSTLES DIIFNSQENL PPLGMKIKKN ANYIELPDLS GDVHPMVQSF AKKVSDVLKQ PQNEAEISKM MPHLKKWQTT LTKLWKQYMN AN
Source	Mammalian cell
Target Names	C03F11.2
Protein Names	Recommended name: Uncharacterized protein C03F11.2
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.