



# Recombinant *Saccharomyces cerevisiae* Protein SLF1 (SLF1)

<b>Product Code</b>	CSB-YP607699SVG
<b>Abbreviation</b>	SLF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12034
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSQNLNDNP KNTSSAAEDK KKQTSSLKLA PIPTTSPWKS SSPDSNTVIP VEELRDISKT AKPSKNGSGS IKLTSNTKWT PITPSVIISG SKDTNSKSGK NSKNSKTNKK MKKRGKYNNND INKKDFNGQT NSTSEISNVS NLESKPLDAN AKVNIHSSSG ATANGNIKRI TNNNNSTNGR QSRNYQNRNG KTRYNNNSRH SQAANNAISF PNNYQARPEY IPNASHWLNN NSRNSYKQLS YFRQQQYYNN INYQQQLQTP YYYSMPIFK SIESIKNQIE FYFSEENLKT DEFLRSKFKK ANDGFIPMSL IGKFYRMVNL SLGGDPNLIL ASMREVLQHK ETNHLEIALG SIEGAQKNMA DDFNPLENYF IRRENWAEYA MESNFDENDD ETEKYNIEKL LGPNDLDNYS YMGYPNFFPS NENGKKSQSY DQGEISRQFE QNLQIND
<b>Source</b>	Yeast
<b>Target Names</b>	SLF1
<b>Protein Names</b>	Recommended name: Protein SLF1
<b>Expression Region</b>	1-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.