



Recombinant *Saccharomyces cerevisiae* Putative uridine kinase DAS2 (DAS2)

Product Code	CSB-BP607705SVG
Abbreviation	DAS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12084
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDRKAVEEKR IVISIGGGHA TGVGAIALDL QNTFKSLYNS INIRVINLDN MIEGNIKSYN NNDYDFDNIL NLVYEKHAVT SQNDMIQHDY EDPIDLIIVC GCYALYDKRI NEISQLKVFL DSDADKRLIS LIKKKNVGSN EQLAQLITEY MDHLRPEMQQ YIEPTRTFAD LIIPSTNENL GRAVLVDGIV KAIEDTKSQI EGNNTNKNIR PRLWDFEAET MDLEKDRYYD LS
Source	Baculovirus
Target Names	DAS2
Protein Names	Recommended name: Putative uridine kinase DAS2 EC= 2.7.1.48 Alternative name(s): DST1-delta 6-azauracil sensitivity protein 2 Regulator of rDNA transcription 3
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.