



Recombinant *Oryza sativa* subsp. japonica Phytoene dehydrogenase, chloroplastic/chromoplastic (PDS)

Product Code	CSB-BP607737OFG
Abbreviation	PDS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0DUI8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	AGQ LSSFFRNSEQ PTKPLQVVIA GAGLAGLSTA KYLADAGHKP ILLEARDVLG GKIAAWKDED GDWYETGLHI FFGAYPNIQN LFGELGINDR LQWKEHSMIF AMPNKPGEFS RFDFPETLPA PLNGIWAILR NNEMLTWPEK VKFALGLLPA MVGGQAYVEA QDGFTVSEWM KKQGVDPDRVN DEVFIAMSKA LNFINPDELS MQCILIALNR FLQEKHGSKM AFLDGNPPER LCMPIVDHVR SLGGEVRLNS RIQKIELNPD GTVKHFALTD GTQITGDAYV FATPVDILKL LVPQEWKEIS YFKKLEKLVG VPVINVHIWF DRKLKNTYDH LLFSRSSLLS VYADMSVTCK EYYDPNRSML ELVFAPAEW VGRSDTEIIE ATMQELAKLF PDEIAADQSK AKILKYHVVK TPRSVYKTIP DCEPCRPLQR SPIEGFYLAG DYTKQKYLAS MEGAVLSGKL CAQSVVEDYK MLSRRSLKSL QSEVPVAS
Source	Baculovirus
Target Names	PDS
Protein Names	Recommended name: Phytoene dehydrogenase, chloroplastic/chromoplastic EC= 1.14.99.- Alternative name(s): Phytoene desaturase
Expression Region	88-578
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.