



Recombinant *Oryza sativa* subsp. japonica Probable transcription factor FL (FL)

Product Code	CSB-EP608017OFG-B
Abbreviation	FL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JAI1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MDPNDAFSAA HPFRWDLGPP APAPVPPPPP PPPPPPPANV PRELEELVAG YGVRMSTVAR ISELGFTAST LLAMTERELD DMMAALAGLF RWDLLLGERF GLRAALRAER GRLMSLGGRH HGHQSGSTVD GASQEVLSDE HDMAGSGGMG DDDNGRRMVT GKKQAKKGS A ARKGKKARRK KVDDLRLDMQ EDEMDCCEDED GGGGSESTES SAGGGGGERQ REHPFVVTEP GEVARAKKNG LDYLFHLYEQ CRLFLLQVQS MAKLHGHKSP TKVTNQVFRY AKKVGASYIN KPKMRHYVHC YALHCLDEEA SDALRRAYKA PGENVGAWRQ ACYAPLVDIS ARHGFIDAV FAAHPRLAIW YVPTLRQLC HQRSSHAAA AAALPPPLF
Source	E.coli
Target Names	FL
Protein Names	Recommended name: Probable transcription factor FL Short name= RFL
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.