



Recombinant *Oryza sativa* subsp. *japonica* Probable esterase PIR7A (PIR7A)

Product Code	CSB-YP608021OFG
Abbreviation	PIR7A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JG98
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEDGGKHFVF VHGLGHGAWC WYRVVAALRA AGHRATALDM AAAGAHPARA DEVGSLEEYS RPLLDAAAA APGERLVLVG HSLGGLSLAL AMERFPDKVA AAVFLAACMP AAGKHMGITL EEFMRRIKPD FFMSDKTIVL NTNQEPRTAV LLGPKLLAEK LYNRSPPEDL TLATMLVRPG TNYIDDPIMK DETLLEGGNY GSVKRVFLVA MDDASSDEEM QRWTIDLSPG VEVEELAGAD HMAMCSKPRE LCDLLLRIAA KYD
Source	Yeast
Target Names	PIR7A
Protein Names	Recommended name: Probable esterase PIR7A EC= 3.1.-.-
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.