



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 11 (Os01g0834700, LOC_Os01g61830)

Product Code	CSB-BP608023OFG
Abbreviation	Os01g0834700, LOC_Os01g61830
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JHZ2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MPPKKAAPSK ADLAKKQKVV EDKTFGLKNK NKSKNVQKYV QSLHQAVQPK PDPTKTAACK KKEEEKAREK ELNDFKVVAV SQPKVPVGVDPK SIVCEFFK VGQCQKGFKC KFSHDLNVQR KGEKIDIYTD KRDAETMEDW DQETLEKVVA SKGAEQQNK PTDIVCKYFL DAVEKKQYGW FWVCPNGGKD CHYRHALPPG YVLKSQMKAL LEEESEKIAI EDEIEDQRKK VKTTTPMTTD LFMEWKRKKA EEREAGLAAL RAERAKNDRM SGRELFMADS SVFVDDAEAY DVYERQEESE ANEEPSNKNQ DEGSSSTSN GKEVEESDDE DINIDDDLDI DELNELEASL SRTSIQIREP GEGTSS
Source	Baculovirus
Target Names	Os01g0834700
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 11 Short name= OsC3H11
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.