



# Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 3 (Os01g0513100, LOC\_Os01g32964)

<b>Product Code</b>	CSB-BP608027OFG
<b>Abbreviation</b>	Os01g0513100, LOC_Os01g32964
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0JMD4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSPSPSSEA AAAHHHHHQR RQHAGAAGGS GLVPLAALIK EEARAERPMG SGSRICARDE EDGGGGGGAE GGRRWRRPLL RYGCAAQSKK GEDFFLLRTD CARPSTSSSS SSSLASSPPH TFAVFAVLDTG HNGNAAAIYT RDNLLNHVLS AMPRGLSREE WLHALPRALV AGFVKTDKEF QHKGQTSGTT ATFVIIDGWT ITVASVGDSR CILDAQGGAV SLLTVDHRLE ENVEERERV ASGGEVGRSL VVGGAIEIGPL RCWPGGLCLS RSIGDIDVGE FIVVPYVKQ VKLSNAGGRL IASDGIWDA LSSEAAAKCC RGLPAELAAK QVVKEALRTR GLKDDTTCIV VDMIPPDQTI RHPSPPKIN KLKSLIFRKK TKDHPNKLTK QLSAAGMVEE LFEEGSAML S ERLGNDSSGR RTSSSLFTCA ICQVDLEPSE GISVHAGSIF SSSSSKPWEG PFLCSDCRDK KDAMEGKRPS GVKVL
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os01g0513100
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 3 Short name= OsPP2C03 EC= 3.1.3.16
<b>Expression Region</b>	1-485
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.