



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 3

(Os01g0513100, LOC_Os01g32964)

Product Code	CSB-EP608027OFG-B
Abbreviation	Os01g0513100, LOC_Os01g32964
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JMD4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MSSPSPSSEA AAAHHHHHQR RQHAGAAGGS GLVPLAALIK EEARAERPMG SGSRICARDE EDGGGGGGAE GGRRWRRPLL RYGCAAQSKK GEDFFLLRTD CARPSTSSSS SSSLASSPPH TFAVFAVLDTG HNGNAAAIYT RDNLLNHVLS AMPRGLSREE WLHALPRALV AGFVKTDKEF QHKGQTSGTT ATFVIIDGWT ITVASVGDSR CILDAQGGAV SLLTVDHRLE ENVEERERV ASGGEVGRSL VVGGAIEIGPL RCWPGGLCLS RSIGDIDVGE FIVPVPYVKQ VKLSNAGGRL IASDGIWDA LSSEAAAKCC RGLPAELAAK QVVKEALRTR GLKDDTTCIV VDMIPPDQTI RHPSPPKIN KLKSLIFRKK TKDHPNKLTK QLSAAGMVEE LFEEGSAML ERLGNDSSGR RTSSSLFTCA ICQVDLEPSE GISVHAGSIF SSSSSKPWEG PFLCSDCRDK KDAMEGKRPS GVKVL
Source	E.coli
Target Names	Os01g0513100
Protein Names	Recommended name: Probable protein phosphatase 2C 3 Short name= OsPP2C03 EC= 3.1.3.16
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.