



# Recombinant *Oryza sativa* subsp. japonica Salt stress root protein RS1 (Os01g0233000, LOC\_Os01g13210)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP608028OFG   |
| <b>Abbreviation</b>      | Os01g0233000, LOC_Os01g13210  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q0JPA6  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. japonica (Rice)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MTSVWTKTKVL TGLNKLFDKD GKKA AAAEFL KSFNKEEIGK EIDDKKTELE PKVVEVVES PPEIKALLKD KKTASKIKKN GPAVTKFLEE LAKIDFPGAK PVSDAVAKSG TTPLSPAIAF ILEKVAPFVP KEEPKPEPEA EAAAETTSRE VAVEEEKKEE EAAPAEPAAA AAEEAAPSTE VVEEKKEEEK PAEAAAPAAE PEKQ   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | Os01g0233000  |
| <b>Protein Names</b>     | Recommended name: Salt stress root protein RS1  |
| <b>Expression Region</b> | 1-204   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |