



Recombinant Danio rerio Protein TSSC4 (tssc4)

Product Code	CSB-BP608096DIL
Abbreviation	tssc4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P4A6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSDRKRKGND APNSLSNRDT VELPDDL SLS DSDSDEPLES FGRKIEDLSS SSEEEGETRQ QVQSVLQENP SFRLTGGSSS FCDRSRDIFA QLDSA AKLTS KDLGEDNILD GTLARPAPPS PPHAVQNKCG VGQELTKKQP PGKKLPDYLA HPERWTHYSL EDIDETSDKK NSQVAHEYIR ELQDSKRSQK AALETFTPAF NQDHGSSNEN KIVFAKPKSK DQSGSRADHA KKEEVGLQHL DDRVEADEGE ALQQSSSWPN KEKKRKWGV TKEEDDTVAPS AVFTSSKRVN RKNFRKTPDD NDGDKE
Source	Baculovirus
Target Names	tssc4
Protein Names	Recommended name: Protein TSSC4
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.