



Recombinant Escherichia coli O6:K15:H31 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase (metE)

Product Code	CSB-EP608364EGY-B
Abbreviation	metE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TAP6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTILNHTLGF PRVGLRRELK KAQESYWAGN STREELLA VG RELRARHWDQ QKQAGIDLLP VGDFAWYDHV LTTSLLLGNV PARHQNKDGS VDIDTLFRIG RGRAPTGEPA AAAEMTKWFN TNYHYMVPEF VKGQQFKLTW TQLLEEVD EALALGHKVKPV LLGPVTYLWL GKVKGQFDR LSLNDILPV YQQVLAELAK RGIEWVQIDE PALVLELPQA WLDAYKPAYD ALQGQVKLLL TTYFEGVTPN LDTITALPVQ GLHVDLVHGK DDVAELHKRL PSDWLLSAGL INGRNVWRAD LTEKYAQIKD IVGKRDLWVA SSCSLLHSPI DLSVETRLDA EVKSWFAL QKCHELALLR DVLNSGDTAA LAEWSAPIQA RRHSTRVHNP AVEKRLAAIT AQDSQR TNVY EVRAEAQRAR FKLPAWPTTT IGSFPQTTEI RTLRLDFKKG NLDANNYRTG IAEHIKQAIV EQERLGLDVL VHGEAERNDM VEYFGEHLDG FVFTQNGWVQ SYGSRVCKPP IVIGDVS RPA PITVEWAKYA QSLTDKPVKG MLTGPVTILC WSPREDVSR ETIAKQIALA LRDEVADLEA AGIGIIQIDE PALREGLPLR RSDWDAYLQW GVEAFRINAA VAKDDTQIHT HMCYCEFNDI MDSIAALDAD VITIETSRSD MELLESEFEF DYPNEIGPGV YDIHSPNVPS VEWIEALLKK AAKRIPAERL WVNPD CGLKT RGWPETRAAL ANMVQAAQNL RRG
Source	E.coli
Target Names	metE
Protein Names	Recommended name: 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase EC= 2.1.1.14 Alternative name(s): Cobalamin-independent methionine synthase Methionine synthase, vitamin-B12 independent isozyme
Expression Region	1-753
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.