



# Recombinant Escherichia coli O6:K15:H31 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase (metE)

<b>Product Code</b>	CSB-MP608364EGY
<b>Abbreviation</b>	metE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TAP6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTILNHTLGF PRVGLRRELK KAQESYWAGN STREELLA VG RELRARHWDQ QKQAGIDLLP VGDFAWYDHV LTTSLLLGNV PARHQNKDGS VDIDTLFRIG RGRAPTGEPA AAAEMTKWFN TNYHYMVPEF VKGQQFKLTW TQLLEEVD EALALGHKVKPV LLGPVTYLWL GKVKGEQFDR LSLNDILPV YQQVLAELAK RGIEWVQIDE PALVLELPQA WLDAYKPAYD ALQGQVKLLL TTYFEGVTPN LDTITALPVQ GLHVDLVHGK DDVAELHKRL PSDWLLSAGL INGRNVWRAD LTEKYAQIKD IVGKRDLWVA SSCSLLHSPI DLSVETRLDA EVKSWFAFAL QKCHELALLR DVLNSGDTAA LAEWSAPIQA RRHSTRVHNP AVEKRLAAIT AQDSQRTNVY EVRAEAQRAR FKLPAPWPTTT IGSFPQTTEI RTLRLDFKKG NLDANNYRTG IAEHIKQAIV EQERLGLDVL VHGEAERNDM VEYFGEHLDG FVFTQNGWVQ SYGSRVCVPP IVIGDVSRPA PITVEWAKYA QSLTDKPVKG MLTGPVTILC WSFPREDVSR ETIAKQIALA LRDEVADLEA AGIGIIQIDE PALREGLPLR RSDWDAYLQW GVEAFRINAA VAKDDTQIHT HMCYCEFNDI MDSIAALDAD VITIETSRSD MELLESEFEF DYPNEIGPGV YDIHSPNVPS VEWIEALLKK AAKRIPAERL WVNPDCLGKT RGWPETRAAL ANMVQAAQNL RRG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	metE
<b>Protein Names</b>	Recommended name: 5-methyltetrahydropteroyltriglutamate--homocysteine methyltransferase EC= 2.1.1.14 Alternative name(s): Cobalamin-independent methionine synthase Methionine synthase, vitamin-B12 independent isozyme
<b>Expression Region</b>	1-753
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.



<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.