



Recombinant Escherichia coli O6:K15:H31 Ribonuclease PH (rph)

Product Code	CSB-BP608369EGY
Abbreviation	rph
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TBG6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRPAGRSNNQ VRPVTLTRNY TKHAEGSVLV EFGDTKVLCT ASIEEGVPRF LKGQQGWIT AEYGMLPRST HTRNAREAAK GKQGGRTMEI QRLIARALRA AVDLKALGEF TITLDCDVLQ ADGGTRTASI TGACVALADA LQQLVENGKL KTNPMKGMVA AVSVGIVNGE AVCDLEYVED SAAETDMNVV MTEDGRIIEV QGTAEGEPFT HEELLTLLAL ARGGIESIVV TQKAALAN
Source	Baculovirus
Target Names	rph
Protein Names	Recommended name: Ribonuclease PH Short name= RNase PH EC= 2.7.7.56 Alternative name(s): tRNA nucleotidyltransferase
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.