



Recombinant Escherichia coli O6:K15:H31 30S ribosomal protein S8 (rpsH)

Product Code	CSB-MP608379EGY
Abbreviation	rpsH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCF5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSMQDPIADM LTRIRNGQAA NKAAVTMPSS KLKVAIANVL KEEGFIEDFK VEGDTKPELE LTLKYFQGKA VVESIQRVSR PGLRIYKRKD ELPKVMAGLG IAVVSTSKGV MTDRAARQAG LGGEIICYVA
Source	Mammalian cell
Target Names	rpsH
Protein Names	Recommended name: 30S ribosomal protein S8
Expression Region	1-130
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.