



Recombinant Escherichia coli O6:K15:H31 GTPase Der (der)

Product Code	CSB-EP608401EGY
Abbreviation	der
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TEX4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MVPVVALVGR PNVGKSTLFN RLTRTRDALV ADFPGLTRDR KYGRAEIEGR EFICIDTGGI DGTEDGVETR MAEQSLLAIE EADVVLFMVD ARAGLMPADE AIAKHLRSRE KPTFLVANKT DGLDPDQAVV DFYALGLGEI YPIAASHGRG VLSLLEHVLL PWMEDLAPQE EVDEDAEYWA QFEAEENGEE EEEDDFDPQS LPIKLAIVGR PNVGKSTLTN RILGEERVVV YDMPGTTRDS IYIPMERDGR EYVLIDTAGV RKRKGITDAV EKFSVIKTLQ AIEDANVVML VIDAREGISD QDLSLLGFIL NSGRSLVIVV NKWDGLSQEV KEQVKETLDF RLGFDIFARV HFISALHGSG VGNLFESVRE AYDSSTRRVG TSMLTRIMTM AVEDHQPLV RGRRVKLKYA HAGGYNPIV VIHGNQVKDL PDSYKRYLMN YFRKSLDVMG SPIRIQFKEG ENPYANKRNT LTPTQMRKRK RLMKHKKSK
Source	E.coli
Target Names	der
Protein Names	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.