



Recombinant Escherichia coli O6:K15:H31 GDP-mannose pyrophosphatase nudK (nudK)

Product Code	CSB-YP608403EGY
Abbreviation	nudK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TF11
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTQQITLVKD KILSDNYFTL HNITYDLTRK DGEVIRHKRE VYDRGNGATI LLYNAKKKT V LIRQFRVAT WVNGNESGQL IETCAGLLDN DEPEVCIRKE AIEETGYEVG EVRKLFE LYM SPGGVTELIH FFIAEYSDSQ RANAGGGVED EDIEVLEL PF SQALEM IKTG EIRDGKT VLL LNYLQMSHLM D
Source	Yeast
Target Names	nudK
Protein Names	Recommended name: GDP-mannose pyrophosphatase nudK EC= 3.6.1.-
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.