



Recombinant Escherichia coli O6:K15:H31 Probable malate:quinone oxidoreductase (mqo)

Product Code	CSB-BP608410EGY
Abbreviation	mqo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TFN1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKKVTAMLFS MAVGLNAVSM AAKAKASEEQ ETDVLLIGGG IMSATLGTYL RELEPEWSMT MVERLEGVAQ ESSNGWNNAG TGHSALMELN YTPQNADGSI SIEKAVAIN E AFQISRQFWA HQVERGVLRT PRSFINTVPH MSFVWGEDNV NFLRARYAAL QQSSLFRGMR YSEDHAQIKE WAPLVMGRD PQQKVAATR EIGTDVNYGE ITRQLIASLQ KKS NFLSLQLS SEVRALKRND DNSWTVTVAD LKNGTAQNIR AKFVFIGAGG AALKLLQESG IPEAKDYAGF PVGGQFLVSE NPDVVNHHLA KVYGKASVGA PPMSVPHIDT RVLDGKRVVL FGPFATFSTK FLKNGSLWDL MSSTTTSNVM PMMHVGLDNF DLVKYLVSQV MLSEEDRFEA LKEYYPQAKK EDWRLWQAGQ RVQIIKRDAD KGGVLR LGTE VVSDQQGTIA ALLGASPGAS TAAPIMLNLL EKVFGDRVSS PQWQATLKAI VPSYGRKLNG DVAATERELQ YTSEVLGLKY DKPQAADSTP KPQLKPQLVQ KEVADIAL
Source	Baculovirus
Target Names	mqo
Protein Names	Recommended name: Probable malate:quinone oxidoreductase EC= 1.1.5.4 Alternative name(s): MQO Malate dehydrogenase [quinone]
Expression Region	1-548
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.