



# Recombinant Escherichia coli O6:K15:H31 Galactose/methyl galactoside import ATP-binding protein MglA (mgIA)

<b>Product Code</b>	CSB-YP608411EGY
<b>Abbreviation</b>	mgIA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TFU2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVSSTTPSSG EYLLEMSGIN KSFPGVKALD NVNLKVRPHS IHALMGENGA GKSTLLKCLF GIYQKDSGTI LFQGKEIDFH SAKEALENGI SMVHQELNLV LQRSVMDNMW LGRYPTKGMF VDQDKMYRET KAIFDELID IDPRARVGT SVSQMQMIEI AKAFSYNAKI VIMDEPTSSL TEKEVNHLFT IIRKLKERGC GIVYISHKME EIFQLCDEV T VLRDGQWIAT EPLAGLTMDK IIAMMVGRSL NQRFPDKENK PGEVILEVRN LTSLRQPSIR DVSFDLHKGE ILGIAGLVGA KRTDIVETLF GIREKSAGTI TLHGKKNH NANEAINHGF ALVTEERRST GIYAYLDIGF NSLISNIRNY KNKVGLLDNS RMKSDTQWVI DSMRVKTPGH RTQIGSLSGG NQQKVIIGRW LLTQPEILML DEPTRGIDVG AKFEIYQLIA ELAKKGGKII IISSEPELL GITDRILVMS NGLVSGIVDT KTTTQSEILR LASLHL
<b>Source</b>	Yeast
<b>Target Names</b>	mgIA
<b>Protein Names</b>	Recommended name: Galactose/methyl galactoside import ATP-binding protein MglA EC= 3.6.3.17
<b>Expression Region</b>	1-506
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.