



Recombinant Escherichia coli O6:K15:H31 Quinate/shikimate dehydrogenase (ydiB)

Product Code	CSB-YP608420EGY
Abbreviation	ydiB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0THD7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDVTAKYELI GLMAYPIRHS LSPMQNKAL EKAGLPFTYM AFEVDNDSFP AAIEGLKALK MRGTGVSM PN KQLACEYVDE LTPAAKLVGA INTIVNDDGY LRGYNTDGTG HIRAIKESGF DIKGKTMVLL GAGGASTAIG AQGAIEGLKE IKLFNRRDEF FDKALAFQR VNENTDCVVT VTDLADQQAF AEALASADIL TNGTKVGMKP LENESLVNDI SLLHPGLLVT ECVYNPHMTK LLQQAQQAGC KTIDGYGMLL WQGAEQFTLW TGKDFPLEYV KQVMGFGA
Source	Yeast
Target Names	ydiB
Protein Names	Recommended name: Quinate/shikimate dehydrogenase EC= 1.1.1.282 Alternative name(s): NAD-dependent shikimate 5-dehydrogenase 2
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.