



Recombinant Escherichia coli O6:K15:H31 Peptidase T (pepT)

Product Code	CSB-BP608428EGY
Abbreviation	pepT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TIU6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDKLLERFLN YVSLDTQSKA GVRQVPSTEG QWKLHLLKE QLEEMGLINV TLSEKGTLM TLPANVPGDI PAIGFISHVD TSPDCSGKNV NPQIVENYRG GDIALGIGDE VLSPVMFPVL HQLLGQTLIT TDGKTLLGAD DQAGIAEIMT ALAVLQQKNI PHGDIRVAFT PDEEVGKGAK HFDVDAFDAR WAYTVDGGGV GELEFENFNA ASVNIKIVGN NVHPGTAKGV MVNALSAAAR IHAEPAPDES PEMTEGYEGF YHLASMKGTV ERADMHYIIR DFDRKQFEAR KRKMMEIAKK VGKGLHPDCY IELVIEDSY NMREKVV EHP HILDIAQQAM RDCDIEPELK PIRGGTDGAQ LSFMGLPCPN LFTGGYNYHG KHEFVTLEGM EKAVQVIVRI AELTAQRK
Source	Baculovirus
Target Names	pepT
Protein Names	Recommended name: Peptidase T EC= 3.4.11.4 Alternative name(s): Aminotripeptidase Short name= Tripeptidase Tripeptide aminopeptidase
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.