



# Recombinant Escherichia coli O6:K15:H31 Peptidase T (pepT)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP608428EGY                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Abbreviation</b>      | pepT                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                     |
| <b>Uniprot No.</b>       | Q0TIU6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen Species</b> | Escherichia coli O6:K15:H31 (strain 536 / UPEC)                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Sequence</b>          | MDKLLERFLN YVSLDTQSKA GVRQVPSTEG QWKLHLLKE QLEEMGLINV<br>TLSEKGTLM TLPANVPGDI PAIGFISHVD TSPDCSGKNV NPQIVENYRG<br>GDIALGIGDE VLSPVMFPVL HQLLGQTLIT TDGKTLLGAD DQAGIAEIMT<br>ALAVLQQKNI PHGDIRVAFT PDEEVGKGAK HFDVDAFDAR WAYTVDGGGV<br>GELEFENFNA ASVNIKIVGN NVHPGTAKGV MVNALS LAAR IHAEPAPDES<br>PEMTEGYEGF YHLASMKGTV ERADMHYIIR DFDRKQFEAR KRKMMEIAKK<br>VGKGLHPDCY IELVIEDSY NMREKVVEHP HILDIAQQAM RDCDIEPELK<br>PIRGGTDGAQ LSFMGLPCPN LFTGGYNYHG KHEFVTLEGM EKAVQVIVRI<br>AELTAQRK |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Target Names</b>      | pepT                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Names</b>     | Recommended name: Peptidase T EC= 3.4.11.4 Alternative name(s):<br>Aminotripeptidase Short name= Tripeptidase Tripeptide aminopeptidase                                                                                                                                                                                                                                                                                                                                                |
| <b>Expression Region</b> | 1-408                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                            |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                      |



of lyophilized form is 12 months at -20°C/-80°C.