



Recombinant Probable small nuclear ribonucleoprotein Sm D1 (snr-3)

Product Code	CSB-EP608718CXY
Abbreviation	snr-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10013
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MKLVRFMLKL SHETVNIELK NGTQVSGTIM GVDVAMNTHL RAVSMTVKNK EPVKLDTL SI RGNNIRYIIL PDPLALDTLL IDDEPRKKAR AARAGASRGR GRGGMRRGGRG GRGRGRGGPR GGGPRR
Source	E.coli
Target Names	snr-3
Protein Names	Recommended name: Probable small nuclear ribonucleoprotein Sm D1 Short name= Sm-D1 Alternative name(s): snRNP core protein D1
Expression Region	1-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.