



Recombinant Saposin-like protein 11 (spp-11)

Product Code	CSB-EP608720CXY-B
Abbreviation	spp-11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10018
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	FV KPQYNVTGQI DSALQRFFGI TLPSLKIPDL LNPDKKRNP SVGQLKKT SF PLCNVNL PPI FFTISLFR IK LPNLIPTALP VIKLPTIKIP NILPTLPTIK VPTIKIPDII PITLPTIKIP EVVPTNLPTV EIPHFIPKTL PTVKIPNIIP TNFPTIETPD IIPKILPTIK IPEIPLTLP TVKIPDIIPI TLPTIKIPEI VPTKLPTVEV PDTIPKTLPT TKIPDIVPIT SPTVKIPQII PTIKIPDIIP KNLSTLGPIK LPTIKLPTLP HILPTKSPKP TPSHKGNMVC DICEKVIGVL TTRLLEIIQK FRVEADKFLT KLCTSLTSNP KTLTVGTM CV MFKGNIMDTI FKGFDGLKKN LEPVSFCKHV PFCK
Source	E.coli
Target Names	spp-11
Protein Names	Recommended name: Saposin-like protein 11
Expression Region	19-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.