



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1D4.01 (SPAC1D4.01, SPAC1F3.11)

Product Code	CSB-EP608746SXV-B
Abbreviation	SPAC1D4.01, SPAC1F3.11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10148
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNSIHIKKKS NRSFRRRKVF GNEKEFDLEE LDDNDIRLRQ ALEATKRRKI RNSIIGINAE KLLNQETKKE KQLNTANEPH EANDQTSQAQS SKLIEAQLPT VEDRFAKQTN EVDINTHLLN FVEKCLKQER LAQNYSENGE TNALNTKNES TVQNIKNSLH PNEHSFIRDA AALGAIREVD LGIISTDVDN LKNGRKRQKK RARMKEKLDLS KALRTSEDA RDEFIEKMLK PISQDEESKG IYRRFRVYKD GTQD
Source	E.coli
Target Names	tls1
Protein Names	Recommended name: Uncharacterized protein C1D4.01
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.