



Recombinant Schizosaccharomyces pombe Meiotically up-regulated gene 113 protein (mug113)

Product Code	CSB-BP608752SXV
Abbreviation	mug113
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10180
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKMPNTLHSV SISTDFPSN VESEAADVNE YSSTSKFETL GNQSSTNLGF SNSLQDKWDC WKRDLSIWW SLSRKATRFY RWLSRSLKWR PVTYETVYPQ TYETQMSEPR EVTTIVKDLN NGDSFVLNVT EPVDPEFLRA NIPPVHRKHL PPRLSLVASP GTIPTRGVVV LEDWINPLLS ERCKLLLQSE LCNQDSYDCP GYICVYRFEQ KNNPSVVLGD STLIQISRVS DLQRHLREYP KNCAFSRSVL EIFPDPGKTS KPCQVSFKVE RLVHKELNEY LSWMQPFTCD SCGSLHENWL QIDSKQWDQI RGVILRWVEY SRVIYA
Source	Baculovirus
Target Names	mug113
Protein Names	Recommended name: Meiotically up-regulated gene 113 protein
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.