



Recombinant Schizosaccharomyces pombe Guanine nucleotide-binding protein subunit beta- like protein (rkp1)

Product Code	CSB-EP608769SXV-B
Abbreviation	rkp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10281
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPEQLVLRAT LEGHSGWVTS LSTAPENPDI LLSGSRDKSI ILWNLVRDDV NYGVAQRRLT GHSHFVSDCA LSFDSHYALS ASWDKTIRLW DLEKGECTHQ FVGHTSDVLS VSISPDNRQV VSGSRDKTIK IWNIIGNCKY TITDGGHSDW VSCVRFSPNP DNLTFFVSAGW DKAVKVDWLE TFLSRTSHYG HTGYVSAVTI SPDGSLCASG GRDGTMLMLWD LNESTHLYSL EAKANINALV FSPNRYWLCA ATGSSIRIFD LETQEKVDEL TVDFVGVGKK SSEPECISLT WSPDGQTLFS GWTDNLIRVW QVTK
Source	E.coli
Target Names	rkp1
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit beta-like protein Alternative name(s): Receptor of activated protein kinase C
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.